



Single Limb, Bain, CPRAM Anesthesia Systems

Unless denoted LF (latex free), Breathing bags incorporated into anesthesia circuits contain natural latex rubber, that may possibly result in an allergic reaction.

Cat. No.	Bag (Liter)	Bleed Valve	Style	Notes	Elbow	Case Qty
Fresh Gas and Transport Systems						
29005	none		Fresh Gas Set-Up		Fresh Gas	10
1456250	2	Y	Transport	Washington T	none	20
1462550	0.5	Y	Transport	Washington T	Gas Sampling	20
Modified Jackson-Rees/Mapleson D Circuits						
14005	0.5		12" Modified Jackson-Rees	Replaced by Cat. No. 1400550	Fresh Gas	20
1400150	1	Y	12" Modified Jackson-Rees		Fresh Gas	20
1400250	2	Y	12" Modified Jackson-Rees		Fresh Gas	20
1400550	0.5	Y	12" Modified Jackson-Rees		Fresh Gas	20
354001	1		12" Modified Jackson-Rees		Fresh Gas	20
354006	0.5	Y	12" Modified Jackson-Rees		Fresh Gas	20
354007	1	Y	12" Modified Jackson-Rees		Fresh Gas	20
7145	3	Y	12" Modified Jackson-Rees		Fresh Gas	20
14101	1		40" Mapleson D	Replaced by Cat. No. 13201-951	Fresh Gas	20
14102	2		40" Mapleson D	Replaced by Cat. No. 13202-951	Fresh Gas	20
354103	3		40" Mapleson D		Fresh Gas	20
354109	3	Y	40" Mapleson D	Replaced by Cat. No. 354103	Fresh Gas	20
Bain Co-axial Partial Rebreathing Circuits						
BD2491	none		60" Co-axial		Gas Sampling	10
BD2494	1		60" Co-axial		Gas Sampling	10
BD2498	3		60" Co-axial		Gas Sampling	10
BD7145	1		33" Co-axial		Gas Sampling	10
BD02492	none		72" Co-axial		Gas Sampling	10
CPRAM Co-axial (Controlled Partial Rebreathing Anesthesia Method)						
13005	0.5		6"	Replaced by Cat. No. 1400550	plain	20
13200	none		40" Co-axial		plain	20
13200-951	none		40" Co-axial		Gas Sampling	20
13201	1		40" Co-axial		plain	20
13201-951	1		40" Co-axial		Gas Sampling	20
13201-952	1		40" Co-axial		GS w/tip	20
13201-BQM	1LF		40" Co-axial		Gas Sampling	20
13202	2		40" Co-axial	Replaced by Cat. No. 13202-951	plain	20
13202-665	2		40" Co-axial	Continuous Sampling Line	plain	20
13202-951	2		40" Co-axial		Gas Sampling	20
13205	0.5		40" Co-axial	Replaced by Cat. No. 13201	plain	20
13300	none		72" Co-axial	Replaced by Cat. No. 13300-951	plain	20
13300-944	none		72" Co-axial	Remote Sampling Connector	plain	20
13300-951	none		72" Co-axial		Gas Sampling	20
13301-934	1LF		72" Co-axial	Remote Sampling Connector	plain	20
13301-951	1		72" Co-axial		Gas Sampling	20
13302	2		72" Co-axial	Replaced by Cat. No. 13202-951	plain	20
13302-951	2		72" Co-axial		Gas Sampling	20
13303-951	3		72" Co-axial		Gas Sampling	20
13400	none		60" Co-axial		plain	20
13400-951	none		60" Co-axial		Gas Sampling	20
13401	1		60" Co-axial	Replaced by Cat. No. 13401-951	plain	20
13401-951	1		60" Co-axial		Gas Sampling	20
13402	2		60" Co-axial	Replaced by Cat. No. 13402-951	plain	20
13402-951	2		60" Co-axial		Gas Sampling	20
13403-951	3		60" Co-axial		Gas Sampling	20
13403-540	3		120" Co-axial		plain	10
13900	n/a		RAD Valve (Respiratory Assist Device)		n/a	1