



## RÜSCH EASYTUBE® DOUBLE LUMEN TUBE

### DESIGNED FOR DIFFICULT TRACHEAL OR ESOPHAGEAL INTUBATION

The Rüsch EasyTube is a dual lumen tube designed for difficult or emergency airway intubation and ventilation. The Rüsch EasyTube creates a viable airway whether it is placed in the trachea or in the esophagus. When placed in the esophagus, the Rüsch EasyTube allows the use of fiber-optic devices, or passage of a suction catheter or tube exchanger. The risk of tracheal trauma is reduced to a minimum due to the single lumen at the distal tip. The complete Rüsch EasyTube set includes tube, syringes pre-filled with air, and an extra length suction catheter. Packaged in a convenient roll-up kit.



### ORDERING INFORMATION

PRODUCT #	SIZE	VERSION	CASE QTY
180041	41 Fr	Large	4
180028	28 Fr	Small	4

When ordering, please always indicate: Product # and French Size



### AT A GLANCE

#### FEATURES

- Single lumen at distal tip reduces potential for trauma
- High-volume low-pressure PVC cuff at distal end allows for minimal wall pressure
- Soft, latex-free cuff to seal the pharynx
- Allows the use of fiber-optic devices or passage of suction catheters and tube exchangers through the second lumen when esophageal intubation is achieved
- Special placement rings clearly indicate a visible depth marking that can be used for either tracheal or esophageal intubation
- Large volume syringe, pre-filled with air for the pharyngeal cuff, and small volume syringe, pre-filled with air for the distal cuff, making intubation quick and easy

To place an order contact customer service at  
**1-866-246-6990**

EasyTube, Rusch and the Rusch logo are registered trademarks of Willy Rusch GmbH. US Patent No. 6,923,176. ©2006 Teleflex Incorporated. All rights reserved. Distributed by Teleflex Medical, a division of Teleflex Incorporated.

L0418-0601