



NEB-U-MASK® SYSTEMS
MULTIPLE THERAPIES, ONE SOLUTION

During an acute asthma attack, a patient may be given oxygen or heliox in combination with pharmacologic therapies. The challenge of this treatment strategy is the lack of a device that allows for the simultaneous delivery of medical gases and aerosolized medication. Low oxygen saturation has been correlated inversely with the rate of hospitalization, which puts caregivers in a difficult situation when faced with the need to remove the oxygen mask to facilitate a neb or MDI treatment.

Introducing the Neb-U-Mask, a combination device that allows for the concurrent delivery of aerosolized medications and high concentrations of oxygen or heliox. The system is composed of an innovative wye design featuring a nebulizer connection and MDI adaptor, a high concentration mask and a 750ml gas reservoir bag. The nebulizer connection features a valved port, which maintains a closed system when the neb is not in use.

Neb-U-Mask is available in both adult and pediatric versions, and is packaged with a Micro Mist® nebulizer and colored gas supply tubing to facilitate ease of use. The device nebulizes at angles up to 90° and features a swivel connector, allowing the system to be used in a manner that is most comfortable for the patient.

The Neb-U-Mask advances patient care and optimizes caregiver efficiency by allowing for an effective treatment strategy in emergency situations. Neb-U-Mask... Multiple therapies, One Solution.

ORDERING INFORMATION

ITEM CODE	DESCRIPTION	CASE QTY
1895	Adult Neb-U-Mask System	20
1896	Pediatric Neb-U-Mask System	20



Latex Free

Rx Only

AT A GLANCE

FEATURES

- FDA approved for heliox and oxygen delivery
- Valved nebulizer port helps maintain a closed system at all times
- 360° swivel connector allows the system to be adjusted to the most comfortable position
- Nebulizer port and MDI adaptor allows treatment to be tailored to the patient's needs
- Color coded gas supply tubing increases safety
- Micro Mist® nebulizer delivers drug at angles up to 90°
- Individually packed 20 to a case

To place an order contact customer service at
1-866-246-6990

Hudson RCI, Micro Mist, and Neb-U-Mask are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. ©2008 Teleflex Incorporated. All rights reserved.

2008-0125 Rev 00